

VascuTherm 3[®]

Thermal Compression Therapy



TREATMENT MODALITIES

Compression Therapy - Alternating/Intermittent programmable compression.

Thermal Therapy - Hot or Cold.

Contrast Therapy - Automatically alternates from hot to cold therapy (20 minutes at 49°F and 10 minutes at 105°F repeating continuously) for pain management. Combined with compression to enhance thermal transfer.

DVT Prophylaxis - Decrease the risk of deep venous thrombosis (DVT). Primarily post surgical.

INDICATIONS

Edema - Reduction of edema associated with soft tissue injuries such as burns, postoperative edema, and ligament sprains.

Lymphedema - Treatment of disorders associated with vascular or lymphatic insufficiency such as Chronic Venous Insufficiency (CVI), venous stasis ulcers, post-mastectomy edema and chronic lymphedema.

Arterial Insufficiency - Treat and assist healing of cutaneous ulceration (wounds), reduce wound healing time, enhance arterial circulation (blood flow), reduce compartmental pressures, reduce edema (swelling), reduce the need for anticoagulant (blood thinning) medications. Aids blood flow back to the heart.

Pain - Cold or Heat therapy for pain management.

A Faster Way to Recover

Features:

- ▲ ICELESS COOLING
- ▲ HEAT THERAPY
- ▲ PRECISE TEMPERATURE CONTROL
- ▲ DESIGNATED COMPRESSION AT THE WOUND SITE
- ▲ DVT PREVENTION THERAPY

Patient Benefits:

- ▲ NO THERMAL TISSUE DAMAGE
- ▲ INCREASED RANGE OF MOTION
- ▲ NO INTERRUPTION IN THERAPY
- ▲ 30% - 40% FASTER RECOVERY
- ▲ REDUCE PAIN WITHOUT THE USE OF MEDICATION



For More Information Please Contact:

TF: 1-800-757-4635
Ph: 780-757-0040
Fx: 780-757-0042

VascuTherm 3[®]

Thermal Compression Therapy



Peak Medical Group Your Recovery Partner

Services:

- ▲ PATIENT INSURANCE INVESTIGATION
- ▲ DELIVERY OF VASCUTHERM 3 PRIOR TO SURGERY
- ▲ PATIENT EDUCATION AND SET UP OF VASCUTHERM 3
- ▲ 24/7 LIVE HELP-LINE
- ▲ POST OPERATION FOLLOW UP
- ▲ PATIENT BILLING AND PICK UP OF DEVICE

Clinic Benefits:

- ▲ EASY INTEGRATION OF VASCUTHERM 3
- ▲ PEAK MEDICAL GROUP HANDLES ALL ASPECTS OF THE RECOVERY PROCESS
- ▲ SURGEON SIGNATURE IS ALL THAT IS REQUIRED

AVAILABLE WRAPS

All wraps available in Sterile and Non-Sterile.

Head	Abdomen	DVT Foot
Face	Cervical	DVT Calf
Breast	Hip	Shoulder
Arm	Foot/Ankle	(Standard, Medium, Large)
Full Arm	Universal U	Knee-Half Leg
Back	AI Foot/Calf	(Standard, Medium, Large)

OPTIONS

- Single Patient Use Wraps
- Travel Bag
- IV Pole Bracket

SPECIFICATIONS

- Weight: 15lbs
- Dimensions: 5.25"W x 8.75"H x 14.25"D
- Hospital Grade Power Cord
- Safety: UL Medical Listing
- FDS 510K

CONTRAINDICATIONS

- Acute Thrombophlebitis
- Presumptive evidence of Congestive Heart Failure
- Pre-existing DVT condition
- Deep Acute Venal Thrombosis (Phlebothrombosis)
- Inflammatory Phlebitis Process
- Episodes of Pulmonary Embolism
- Pulmonary Edema
- Acute Inflammations of the Veins (Thrombophlebitis)
- Decompensated Cardiac Insufficiency
- Arterial Dysregulation
- Erysipelas
- Carcinoma and Carcinoma Metastasis in the affected Extremity
- Decompensated Hypertonia
- Acute Inflammatory Skin Diseases or Infection
- Venous or Arterial Occlusive Disease
- Venous or Lymphatic Return is Undesirable
- Severe Arteriosclerosis or active infection
- Sepsis in the limb
- Cellulitis without appropriate antibiotic coverage
- Immediately following skin grafts in/around treatment sites